



STRATOS®

RETRA™

ENDOSCOPIC CARPAL AND CUBITAL TUNNEL RELEASE SYSTEM

- Disposable System Requiring No Additional Instrumentation
- "Static" Cannula to Prevent Subluxation of Pertinent Anatomy
- Clear Cannula for Enhanced Endoscopic Views
- Compatible with Stryker® 4mm, 30° Arthroscope
- Actuating Blade Cuts Retrograde or Antegrade ●●●



A.M. SURGICAL, INC.™

For more information or to order contact us at
800.437.9653 • info@amsurgical.com
www.amsurgical.com