

# SECURE Screw Fixation System for Hand Fractures



Retrograde • Intramedullary • Headless Screw Fixation

## FEATURES

- Minimally invasive approach reduces OR time
- Early ROM and return to ADL
- Obviates prolonged postoperative immobilization
- Extended screw lengths (35mm to 60mm via 5mm increments) to capture proximal diaphysis without marked countersinking
- Tapered T8 driver minimizes potential for articular cartilage injury
- Differential screw pitch facilitates reduction of axially stable diaphyseal fractures without over-compression



## AVAILABLE LENGTHS

35mm

40mm

45mm

50mm

55mm

60mm

\*Screws are individually sterile packaged

For more information or to order contact us at  
800.437.9653 • [info@amsurgical.com](mailto:info@amsurgical.com)  
[www.amsurgical.com](http://www.amsurgical.com)

## INSTRUMENT TRAY INCLUDES

Tapered T8 Cannulated Driver



2.0 mm Cannulated Drill Bit



Measuring Device



1.1mm K-wire (x2)